

6_Therapy


11 March 2025 14:27

Shock Wave Lithotripsy (SWL)

- **Introduction & Adoption:**
 - SWL was introduced in the 1980s as a treatment for urinary stones.
 - It quickly became the first-line treatment option.
 - In the U.S., approximately **70% of kidney stones** are treated using SWL.
- **Technological Advancements:**
 - SWL has seen several waves of technological improvements.
 - However, the **fundamentals of shock wave generation and delivery remain unchanged.**
 - Changes in lithotriptors have been significant from a **user perspective** but not in the **core pressure pulse mechanism.**
- **Mechanism of Action & Side Effects:**
 - Lithotriptors produce an **acoustic shock wave**, a characteristic pressure pulse.
 - This pulse is responsible for **breaking stones** but can also cause **collateral tissue damage**, sometimes significant.

Mechanisms of Shock Wave Lithotripsy (SWL) and Technological Developments

- **Mechanical Forces in SWL:**
 - Lithotriptors generate a **powerful acoustic field**, causing two primary mechanical forces:
 - **Direct stress** from the high-amplitude shock wave.
 - **Stresses and microjets** from the growth and collapse of **cavitation bubbles.**
 - Despite advancements in research, the **exact mechanisms of shock wave action** are not fully understood.
- **Challenges in Lithotripter Development:**
 - No **completely safe yet effective** lithotripter has been developed.
 - Recent trends in **high-amplitude, tightly focused shock waves** have increased **adverse effects and re-treatment rates.**
 - The goal within the field is to make SWL both **safer and more effective.**
- **Understanding Lithotripter Function:**
 - Urologists must **optimize the use of available machines** by understanding how they function.
 - This chapter introduces **fundamental physical concepts** underlying SWL mechanisms.
- **Key Topics Covered in the Chapter:**
 - **Characteristics of lithotripter shock waves**
 - **Acoustics of SWL** (including an **acoustics primer**)
 - **Acoustic cavitation** and its effects
 - **Physics of clinical lithotriptors**
 - **Shock generation and focusing**
 - **Coupling shock waves to the body**
 - **Focal zone of high acoustic pressure**
 - **Mechanisms of shock wave action:**
 - How **shock waves break stones**
 - How **shock waves cause tissue damage**
 - **Evolution of lithotripter technology**
 - **Future directions in lithotripsy**




Characteristics of a Lithotripter Shock Wave

- **Shock Wave Profile (Figure 38-1A):**
 - A typical shock wave at the focus of a lithotripter is a **short pulse** lasting **~5 μs**.
 - The wave starts with an **instantaneous jump** to a **peak positive pressure of ~40 MPa**.
 - This fast transition, referred to as a “**shock**,” is **faster than can be measured (<5 ns)**.
 - The pressure drops to zero **after ~1 μs**, followed by a **negative pressure phase (~3 μs, -10 MPa peak)**.
 - Unlike the positive phase, the **negative phase lacks a sharp transition**.
 - The entire **5 μs pulse** is often called a **shock wave, shock pulse, or pressure pulse**, though **technically only the leading sharp transition is a shock**.
- **Frequency Spectrum (Figure 38-1B):**
 - The **amplitude spectrum** of the shock pulse reveals **no dominant frequency** or tone.
 - The energy is spread **over a wide frequency range** (100 kHz – 1 MHz), a key characteristic of a **short pulse**.
 - This suggests that lithotriptors **do not break kidney stones via resonance excitation**, unlike how an opera singer shatters a crystal glass.

Variations in Lithotripter Shock Waves

- **General Characteristics of Lithotripter Shock Waves:**
 - Most lithotriptors produce a **similarly shaped shock wave**, consisting of:
 - A **shock front**
 - A **compressive phase**
 - A **tensile tail**
 - Despite the similarity in waveform shape, **amplitude varies significantly** depending on the machine and settings.
 - **Peak positive pressure** typically ranges between **30 MPa and 110 MPa**, while **negative pressure** ranges between **-5 MPa and -15 MPa**.
- **Comparison of Shock Waveforms (Figures 38-2A & 38-2B):**



- **Figure 38-2A:**
 - Shows a comparison between **electrohydraulic** and **electromagnetic lithotripter waveforms**.
 - The **main difference is amplitude**, while the basic waveform structure remains similar.
- **Figure 38-2B:**
 - Displays waveforms at **lower power settings** for both lithotriptors.
 - Again, the **shape is preserved**, but the **amplitudes differ**.
- **Impact of Shock Wave Amplitude & Focal Zone:**
 - Lithotripter **shock waves have a unique form**, characterized by:
 - A **high-amplitude compressive phase** with an **extremely rapid transition** and **short duration**.
 - A **trailing tensile phase**.
 - While all lithotriptors share this basic waveform, there are **considerable differences in amplitude and spatial extent** of the acoustic output.
 - **Performance differences** among lithotriptors are likely influenced by **amplitude and focal zone size**.

What is an Acoustic Wave?

- **Definition & Formation:**
 - An **acoustic wave (sound wave)** is created whenever an **object moves within a fluid** (gas or liquid).
 - As the object moves, it **compresses** the surrounding fluid, forcing molecules closer together.
 - These compressed molecules push against **adjacent molecules**, transferring the compression forward.
 - This **chain reaction** results in a **traveling wave of compression**, known as an **acoustic wave**.
- **Wave Propagation & Medium Dependence:**
 - The **speed of an acoustic wave** (sound speed) is a **material property** of the medium.
 - **Key distinction: Molecules do not travel with the wave**—they only **jostle adjacent molecules** to pass energy forward.
 - A **medium is required** for acoustic waves to propagate.
- **Difference Between Acoustic and Electromagnetic Waves:**
 - **Acoustic Waves** (e.g., sound, seismic, water waves) require a **medium** (fluid or solid) to support **vibrations**.
 - **Electromagnetic Waves** (e.g., light, radio waves, X-rays) do **not** require a medium; they are carried by **photons** and can travel through a **vacuum**.
 - This explains why **sound cannot travel in space**, but **light can**.

Compressive and Tensile Phases of Sound Waves

- **Two Phases of an Acoustic Wave:**
 - A **compressive phase** occurs when molecules are **pushed together** due to the movement of an object within a fluid.
 - A **tensile phase (rarefaction)** occurs when an object moves away, creating a **partial vacuum**, causing neighboring molecules to move in and create a propagating **rarefactional disturbance**.
 - **Both compressive and tensile waves propagate at the same sound speed** in most cases.
- **Wave Generation and Form:**
 - **Typical acoustic sources** (e.g., audio speakers) **vibrate back and forth**, producing alternating **compression and rarefaction waves**.
 - This results in a **waveform** that is often **sinusoidal**, but **most acoustic waves, including lithotripter pulses, are not sinusoidal**.
- **Wave Speed and Amplitude:**
 - For **small-amplitude waves (linear acoustics)**, every point of the waveform moves at the **same speed (c_0)**.
 - In **water and tissue**, the **sound speed is $\sim 1,500$ m/s**.
 - For **large-amplitude (nonlinear) waves**, such as **shock waves**, the **sound speed is slightly modified** by the wave's presence.
- **Measuring Acoustic Waves:**
 - Acoustic waves are typically measured as **pressure pulses over time at a fixed point in space** (e.g., via a microphone).
 - However, **waves also vary in space**, and it is useful to consider their **spatial extent** for analysis.

Reflection and Transmission of Sound Waves

- **Acoustic Wave Interaction with Different Media**
 - When an acoustic wave encounters a **medium with different impedance**, part of the wave is **transmitted** into the new medium, and part is **reflected** back.
 - At **normal incidence** (wave perpendicular to the surface), the amplitudes of the transmitted and reflected waves depend on the **impedance mismatch** between the two media.
- **Transmission Coefficient and Impedance Mismatch (Figure 38-4)**
 - The **intensity transmission coefficient** determines how much energy passes into a new medium.
 - **Water-to-tissue transmission:** Very efficient.
 - **Water-to-kidney stone transmission:** High, with **75–95% of energy** transmitted.
 - **Water-to-air transmission:** Extremely low (**$<0.1\%$**), meaning **99.9% of the energy is reflected**.
- **Implications for Lithotripsy**
 - **Shock wave generators are water-filled** to ensure efficient energy transfer.
 - **Immersion of the patient in water** provides the best coupling of shock waves to the body.
 - **Dry lithotriptors require careful air pocket elimination** between the shock head and the body.
 - **Shock waves do not effectively pass through air-filled organs** like the lungs or bowel segments.
 - **Optimal shock wave path:** The **flank of the patient** (defined by the **ribs, spine, and pelvic bone**) provides the best **pure tissue path to the kidney** for effective treatment.

Focusing and Diffraction of Sound

- **Purpose of Focusing in Lithotripsy**
 - **Shock waves are focused** to concentrate **acoustic energy** on the **kidney stone** while minimizing damage to surrounding tissue.

- Lithotriptors achieve focusing through:
 - **Reflectors**
 - **Acoustic lenses**
 - **Spherically curved sources**
- **Wave Focusing Physics**
 - The **fundamental physics** of wave focusing is similar across different methods.
 - **Ideal focus**: All energy concentrated in an **infinitesimally small region** (theoretically impossible).
 - **Diffraction limitation**: Energy cannot be arbitrarily confined due to **wave propagation constraints**.
- **Focal Zone**
 - **Definition**: The finite region around the focus where **acoustic pressure remains high**.
 - **Theoretical lower bound**: If sound approaches from all angles, diffraction limits the **focal zone** to about **one wavelength**.
 - **Practical lithotripsy systems**: Since sound comes from a **single direction**, the **focal zone ranges from a few millimeters to tens of millimeters** in size.

Focal Zone

- **Definition and Terminology**
 - The **focal zone** of a lithotripter is an **ellipsoidal** region where **high acoustic pressure** is concentrated.
 - Equivalent terms:
 - **Focal region**
 - **Hot spot**
 - **Focal spot**
 - **Focal volume**
 - **Zone of high pressure**
 - The **longest dimension** of the focal zone aligns with the **shock wave axis**.
- **Factors Determining Focal Zone Size**
 - The **length and diameter** of the focal zone depend on:
 - **Source diameter**
 - **Focal length** (distance from source to target)
 - **Frequency content** of the waveform
 - **Design features** of different lithotriptors affect their **focal zone size and acoustic pressure distribution**.
 - Some lithotriptors generate **extreme acoustic pressures in very narrow focal zones**.
- **Focal Zone Estimation**
 - **No explicit formula** for pulsed lithotripsy waveforms.
 - The **focal region of a lithotripter** can be approximated using equations for the **focal region of a sine wave**.
- **Focal Zone Adjustments (Figures 38-6A & 38-6B)**
 - **Increasing source aperture diameter** → **Shorter and narrower focal zone**.
 - **Decreasing focal length (source-to-target distance)** → **Shorter and narrower focal zone**.
 - **Optimal design for a small focal zone**:
 - **Large aperture diameter**
 - **Short focal length**
 - **Limitations**:
 - The **size of the acoustic window in the flank** restricts extreme modifications.
 - The lithotripter must still be able to **target stones deep in the body**.


Sound Waves Are Not Just Pressure Waves

- **Multiple Effects of Sound Wave Propagation:**

- A propagating **sound wave** affects more than just **pressure**—it also influences:
 - Density (ρ):**
 - **Increases** in compressed regions.
 - **Decreases** in rarefied (tensile) regions.
 - Pressure:**
 - **Rises** during compression.
 - **Falls** during rarefaction.
 - Particle velocity:**
 - Individual **fluid particles move** in response to the wave, but they **do not travel with it**.
- **Wave Behavior in a Fluid:**
 - Compression causes **higher local density**, while rarefaction leads to **lower local density**.
 - These variations contribute to the overall **wave propagation** in the medium.

Nonlinear Effects in High-Amplitude Acoustic Waves

- **Wave Speed Variation in High-Amplitude Waves:**
 - In **large-amplitude acoustic waves**, such as **lithotripter shock waves**, the **wave speed is no longer constant**.
 - The speed of the wave **depends on the local compression** of the fluid.
- **Nonlinearity in Shock Waves:**
 - Even in **high-power lithotriptors**, water compression is **less than 5%**, yet nonlinearity still occurs.
 - Two physical effects contribute to this nonlinearity:
 - Higher Pressure Increases Local Sound Speed:**
 - In regions of **high pressure**, the **local sound speed increases** beyond the normal value.
 - Convective Effects in High-Pressure Regions:**
 - Molecules in **high-pressure regions** move with **higher particle velocity**, shifting forward in the **direction of wave propagation**.
 - In **tissue**, the **first effect (increased local sound speed)** dominates nonlinearity.
- **Wave Distortion and Nonlinear Propagation:**
 - **Nonlinear waves** differ from **linear waves** because **different parts of the wave travel at different speeds** (Equation 38-19).



- **Figure 38-7** illustrates how a **sinusoidal wave deforms** as it propagates under nonlinear conditions.
- Without **absorption**, the **wave can develop an infinite slope**, causing it to **fold over** and become **multivalued**.
- **This is similar to ocean waves fanning up a beach**, but in acoustics, such a waveform is **not physically possible**, as a **single point in space cannot have multiple pressure values simultaneously**.

Hydrophones Measure Shock Waves

- **Lithotripter Acoustic Field Measurement**
 - The main physical property of a lithotripter is its spatial and temporal acoustic pressure distribution.
 - The acoustic field is typically measured in water using a **hydrophone**, which converts pressure into an electrical signal.
- **Characteristics of Lithotripsy Shock Waves**
 - Lithotriptors generate **short, high-amplitude acoustic pulses** with a wide frequency band.
 - These pulses are focused on a **small volume in space**.
- **Hydrophone Requirements**
 - **Wide bandwidth**: 60 kHz to >20 MHz.
 - **Robustness**: Must withstand high pressures from shock waves.
 - **Small active area**: Approximately **0.5 mm**.
- **PVDF Hydrophones**

- First reliable lithotripsy shock wave measurements were done using **PVDF (polyvinylidene fluoride)**, a **piezoelectric plastic**.
- PVDF has:
 - **Very wide bandwidth**
 - **High amplitude acoustic pressure measurement capability**
 - **Manufacturable with a small active region**
- **Types of Hydrophones**
 - **Membrane hydrophones** and **needle hydrophones** have been used.
 - **Membrane hydrophones** provide the **best shock wave measurement**.
- **Limitations of PVDF Hydrophones**
 - **Weak water adhesion:** The tensile phase of the lithotripsy pulse can cause **cavitation** at the PVDF surface.
 - **Consequences of cavitation:**
 1. **Measurement issues:** The hydrophone cannot measure the **tensile phase** accurately. When a bubble forms, the **negative pressure is relieved**, causing the hydrophone to register close to zero.
 2. **Physical damage:** **Collapsing cavitation bubbles** can irreversibly damage the hydrophone.

Fiber Optic Probe Hydrophone (FOPH) for Lithotripsy Shock Waves

- **Introduction of FOPH**
 - A new hydrophone, the **fiber optic probe hydrophone (FOPH)**, has been developed.
 - It is now considered **state-of-the-art** for measuring lithotripsy shock waves.
 - **International measurement standards** recommend using the FOPH.
- **Working Principle**
 - The FOPH consists of a **laser** that injects light into one end of an **optical fiber**.
 - The other end is placed in the **lithotripter field**.
 - It measures the **light reflected from the fiber's end**, which varies depending on the pressure in the fluid.
- **Advantages Over PVDF Hydrophones**
 - **Wide bandwidth** and **high-pressure measurement capability**, similar to PVDF.
 - **Smaller active area: 100 μm** (compared to **500 μm** for PVDF hydrophones).
 - **Stronger water adhesion:**
 - Made of **silica**, which has high adhesion with water.
 - **Less cavitation**, allowing for **better measurement of the tensile phase**.
 - **More resistant to cavitation damage** than PVDF.
- **Limitations of FOPH**
 - **Weak signal output**.
 - **Not suitable for measuring low pressures** (≈ 2 MPa or less).

Measuring Shock Waves in the Body

- **Current Measurement Limitations**
 - **All published focal zone dimensions** come from **in vitro** experiments.
 - **Few studies** have attempted to measure **shock wave pressure in living animals**.
- **Effects of Tissue on Shock Wave Propagation**
 - **Tissue properties** (high attenuation and inhomogeneity) alter shock wave behavior as they travel through the body.
 - **Comparison of in vitro and in vivo waveforms (Figure 38-9):**
 - Data from **Dornier HM3 lithotripter** measured in **water** and in a **pig using a PVDF membrane hydrophone**.
 - **Basic waveform shape remains similar** between in vitro and in vivo conditions.
 - **Key differences in vivo:**
 - **30% decrease in peak positive pressure**.

- **Significantly increased shock rise time (70 ns).**
- **Both changes are consistent with the increased attenuation of tissue.**

Acoustic Cavitation

- **Definition and Mechanism**
 - **Acoustic cavitation** is a mechanical force generated by lithotripter shock waves.
 - It occurs when the **tensile phase (negative pressure) of an acoustic wave** is strong enough to **create cavities (bubbles) in a fluid**.
 - In lithotripsy, the **tensile phase (~10 MPa)** induces **violent cavitation events**.
- **Role of Cavitation in Lithotripsy**
 - **Tissue damage** during **shock wave lithotripsy (SWL)** is partly caused by cavitation.
 - **Stone comminution** (fragmentation) is also **enhanced** by cavitation effects.
- **Bubble Formation and Growth**
 - Cavitation is typically **initiated** at:
 - **Micron-sized motes** in the fluid.
 - **Small gas pockets** trapped on rough surfaces.
 - Several **theories** describe how cavitation progresses once a bubble forms.
- **Bubble Dynamics (Figure 38-10)**
 - A **spherical bubble's radius** changes in response to a **lithotripter shock wave**:
 - 1. Compression phase:** The **positive pressure** of the shock wave **compresses the bubble**.
 - 2. Expansion phase:** The **tensile phase** causes the bubble to grow from **1 μm to ~1 mm** over **150 μs** .
 - 3. Inertial cavitation:**
 - The bubble continues to **expand after the shock wave passes** (**~5 μs**).
 - At this stage, bubble dynamics are **driven by fluid inertia** rather than acoustics.
 - **Gas and vapor diffusion** occurs inside the bubble.
 - 4. Collapse phase:**
 - The **near-vacuum inside the bubble** and the surrounding **ambient pressure (~1 atm)** cause the bubble to collapse.
 - The collapse is **extremely violent**, generating **high heat and compression** that can even produce **light emission**.
 - 5. Post-collapse effects:**
 - **Bubble rebounds** occur.
 - The **gas inside the bubble** slowly **diffuses back into the fluid**.
- **Acoustic Emissions from Cavitation**
 - **Lithotripsy-induced cavitation bubbles** generate **two acoustic emissions**:
 - 1. Upon impact by the compressive shock wave.**
 - 2. Upon collapse, hundreds of microseconds later.**
 - This **"double-bang" acoustic signature** can be used to **detect cavitation events**.
- **Cavitation Measurement Techniques**
 - Multiple methods exist to **measure cavitation**, utilizing the unique **acoustic emissions** it generates.

High-Speed Photography of Bubble Dynamics

- **Observing Cavitation with High-Speed Cameras**
 - High-speed photography can be used **in vitro** to observe **bubble behavior** from **formation to collapse**.
 - In **theory**, this method allows **tracking the entire bubble lifecycle**, but **practical limitations** exist.
- **Challenges in Capturing Full Bubble Dynamics**
 - **Growth Phase:**
 - Requires imaging at **millimeter scales** and **tens of microseconds** time resolution.

- **Collapse Phase:**
 - Bubble shrinks to **less than 1 μm** .
 - Collapse dynamics occur at **nanosecond** timescales.
- **Rebound and Dissolution Phase:**
 - **Remnant bubble ($\sim 10 \mu\text{m}$)** dissolves slowly over **hundreds of milliseconds**.
- **Wide range of spatial and temporal scales** makes it **impossible** to photograph the full process in a single recording.
- **Segmented Approach to Cavitation Study**
 - Due to these limitations, **cavitation is studied in separate phases** rather than as a whole.
- **Depth of Field Limitations**
 - **Cameras have a limited depth of field**, making it difficult to capture bubble behavior **throughout the entire cavitation field**.
 - This **prevents a complete record** of bubble dynamics in **large cavitation volumes**.

Laser Scattering of Single Bubbles

- **Principle of Laser Scattering for Bubble Measurement**
 - A **laser beam** is used to **illuminate** a cavitation bubble.
 - A **photodetector** collects the **scattered light**, which varies based on the **bubble's radius**.
 - For a **spherical bubble**, the amplitude of scattered light follows a **known pattern**, allowing **precise measurement of bubble dynamics**.
- **Advantages of Laser Scattering**
 - Captures **most spatial and temporal scales** of **lithotripsy-induced cavitation bubbles**.
 - Provides **high precision** for measuring **single bubble dynamics**.
- **Limitations of Laser Scattering**
 - **Small sample volume:** Can only measure a **tiny region** of the cavitation field.
 - **Requires unobstructed visual access** at **high magnification**.
 - **Assumes a single spherical bubble:**
 - Theory is **only accurate for isolated, perfectly spherical bubbles**.
 - **Cannot provide quantitative data** for **bubble clouds** or **non-spherical bubbles**, both of which are **common in lithotripsy-induced cavitation**.
 - Measurements in realistic **cavitation fields** are mostly **qualitative**.

Acoustic Detection Can Be Used In Vivo

- **Advantages of Acoustic Detection**
 - **Powerful tool** for characterizing **bubble dynamics in living subjects**.
 - Unlike imaging techniques, it can **monitor cavitation events in real-time**.
- **Modes of Acoustic Detection**
 - Active Cavitation Detection (ACD)**
 - Works similarly to **laser scattering** but uses **acoustic waves**.
 - **One transducer** sends an **acoustic wave** into the cavitation field.
 - **A second transducer** picks up **reflections from bubbles**.
 - Passive Cavitation Detection (PCD)**
 - **No external wave is transmitted**; instead, transducers **listen for natural emissions** from bubbles.
 - Cavitation bubbles produce a characteristic **"double-bang" acoustic signal**.
 - **Dual PCD:**
 - Uses **two receiving transducers**.
 - Enables **coincidence detection**, allowing sampling of a **small, discrete volume** of the cavitation field.
 - **Signal timing and amplitude** are influenced by factors like:
 - **Initial bubble size**.

- **Lithotripter pulse amplitude.**
- **Limitations of Acoustic Detection**
 - **Does not image bubbles** (cannot measure **bubble number or exact size**).
 - Instead, provides **valuable data** for:
 - **Characterizing lithotripter acoustic output.**
 - **Assessing cavitation field dynamics.**

Other Cavitation Measurement Techniques (Figure 38-11)

- **Metal Foil Pitting Method**
 - Cavitation leads to **pitting on metal foils**.
 - **Number and depth of pits** can indicate the **violence of cavitation**.
- **Electromagnetic Probe**
 - Measures the **mechanical force** on a **steel ball** from:
 - **The incident shock wave.**
 - **Cavitation activity.**
- **Sonoluminescence and Sonochemistry**
 - **Extreme pressures and temperatures** inside cavitation bubbles can:
 - **Produce light emissions (sonoluminescence).**
 - **Enhance chemical reaction rates (sonochemistry).**
 - **Production of light and chemical byproducts** can be used to **quantify cavitation activity**.
 - **Interpretation is complex** as these are **secondary measurements** of the cavitation field.

Shock Generation and Focusing

Laser Scattering of Single Bubbles

- **Principle:**
 - A **laser beam** illuminates a **single cavitation bubble**, and a **photodetector** collects the **scattered light**.
 - The **amplitude of scattered light** varies predictably with **bubble radius**, allowing precise measurement.
- **Advantages:**
 - Captures **most spatial and temporal dynamics** of lithotripsy-induced cavitation bubbles.
- **Limitations:**
 - **Small sample volume** restricts large-scale measurement.
 - Requires **unrestricted visual access at high magnification**.
 - **Assumes a single spherical bubble**, making it **ineffective for bubble clouds or non-spherical bubbles**, which are common in lithotripsy.

Acoustic Detection Can Be Used In Vivo

- **Modes of Acoustic Detection:**
 - Active Cavitation Detection (ACD)**
 - **One transducer** sends an **acoustic wave** into the cavitation field.
 - **A second transducer** detects **reflections from bubbles** (analogous to laser scattering).
 - Passive Cavitation Detection (PCD)**
 - **No external wave transmission**—transducers **listen for natural emissions** from cavitation bubbles.
 - **Dual PCD:**
 - Uses **two receiving transducers**.
 - Allows **coincidence detection**, sampling a **small, discrete cavitation volume**.
- **Benefits of Acoustic Detection:**
 - Provides **valuable data** for **shock wave characterization** and **cavitation field dynamics**.
 - Can be used **in vivo**, unlike optical methods.
- **Limitations:**

- Cannot **image bubbles** (does not provide bubble count or exact size).

Electrohydraulic Lithotriptors (EHL)

- **Shock Wave Generation:**
 - **Spark source** creates a **shock wave**.
 - The wave is **focused by an ellipsoidal reflector** (Figure 38-12).
 - Unlike other lithotriptor types, the **shock wave remains a shock wave throughout propagation**.
- **Focusing Mechanics:**
 - The **spark source must be precisely placed at the first focus of the ellipse**.
 - **Misalignment (even by a few mm)** causes:
 - **Loss of focusing efficiency.**
 - **Lengthening and broadening of the focal zone.**
 - **Electrode variability** causes **shot-to-shot variations**:
 - Shock wave amplitude can vary by **up to 50%**.
 - **Focal zone position may shift.**
- **Unique "Two-Pulse" Feature:**
 - EHL insonifies the target with **two waves**:
 - Direct wave:**
 - Travels **directly from the spark** to the target.
 - Arrives **~30 μs before the main wave**.
 - Has **low amplitude due to spherical spreading**.
 - Can **influence cavitation** caused by the main focused wave.
 - Main focused wave:**
 - Reflected and **concentrated at the focal point**.
- **Electrode Wear and Maintenance:**
 - **Electrodes degrade over time** and require replacement.
 - Some manufacturers **extend electrode life** by:
 - **Encapsulation.**
 - **Filling the casing with electrolyte.**
 - **Electrode wear affects acoustic output**, leading to inconsistencies in shock wave generation.

Electromagnetic Lithotriptors (EML)

- **Shock Wave Generation Mechanism**
 - Uses an **electrical coil** near a **metal plate** as the **acoustic source**.
 - When **electrical current pulses through the coil**, the **plate experiences a repulsive force**, generating an **acoustic wave**.
- **Types of Wave Generation (Figure 38-13)**
 - **Flat metal plate** → Generates a **plane wave**, focused by an **acoustic lens** (Figure 38-13A).
 - **Tube-shaped plate** → Generates a **cylindrical wave**, focused by a **parabolic reflector** (Figure 38-13B).
- **Advantages Over Electrohydraulic Lithotriptors (EHL)**
 - **More consistent focusing:**
 - Variation in **measured pressure waves** is **<10%**.
 - **Less shot-to-shot variation** compared to EHL.
 - **No electrodes to replace**, reducing maintenance needs.
- **Acoustic Differences Between EML and EHL**
 - Initial Waveform Shape:**
 - **EHL:** The wave **starts as a shock wave**.
 - **EML:** The wave starts as a **high-intensity ultrasonic wave** (smooth waveform, no discontinuities).
 - **EML shock formation:** The ultrasonic wave **undergoes nonlinear distortion** during propagation, **becoming a shock wave** before reaching the focus.

Trailing Positive Pressure Difference:

- **EML waveforms have a smaller trailing positive pressure** after the negative phase.
- This likely does not impact stone stress but may influence cavitation dynamics.

Piezoelectric Lithotripter (PEL)

- **Shock Wave Generation Mechanism**
 - Uses **piezoelectric crystals** to generate **ultrasonic waves**.
 - When **voltage is applied**, the **crystals deform**, producing an **acoustic wave**.
 - **Crystals are arranged on the inside of a spherical cap**, focusing the acoustic wave at the **center of the sphere's curvature** (Figure 38-14A).
- **Advantages of PEL**
 - **Highly reproducible focus with minimal variation in focal waveforms**.
 - **Consistent shock wave formation**, similar to **electromagnetic lithotriptors (EML)**.
- **Waveform Characteristics**
 - **Starts as an acoustic pulse**, similar to EML.
 - A **shock wave forms before reaching the focus** due to **nonlinear propagation distortion** (Figure 38-14B).
 - **Presence of a trailing coda (ringing effect):**
 - **Piezoelectric crystals continue vibrating for a few cycles** after excitation.
 - This **ringing effect is not present in EHL or EML**.
 - The **coda does not impact stone stress**, but it **may influence cavitation dynamics**.

Coupling of the Shock Source to the Body

- **Importance of Acoustic Impedance Matching**
 - Efficient **acoustic energy transfer** occurs when **acoustic impedances of two media are similar**.
 - **Water/tissue interface** provides **excellent coupling**, theoretically allowing **>99% energy transfer** into the body.
- **Impact of Air Pockets on Energy Transfer**
 - Even a **small air pocket** at the skin surface **drastically reduces energy transfer**.
 - **Air acts as a strong reflector**, preventing shock waves from penetrating the body effectively (Figure 38-4).
- **Significance of Proper Coupling**
 - Ensuring **good contact between the shock wave source and the patient** is **critical** for effective lithotripsy.

Water-Bath and Dry Lithotriptors

Water-Bath Lithotriptors (First-Generation, e.g., Dornier HM3)

- **Design and Coupling:**
 - **Patient is fully immersed in an open water bath**.
 - **Nothing but water** between the **shock source** and the **patient**, ensuring **ideal acoustic coupling**.
- **Challenges:**
 - **Bubbles from the spark gap and cavitation bubbles** can accumulate **against the patient's skin**, interfering with **subsequent shock waves**.
 - To prevent this:
 - The **ellipsoidal reflector** in the **Dornier HM3** is **tilted 14 degrees off vertical**.
 - **Continuous degassing** of the water is performed to **remove bubbles**.

Dry Lithotriptors (Modern Lithotriptors)

- **Design and Coupling:**
 - The **shock wave source** is mounted in a **“therapy head”** filled with **water**.

- A **thin rubber membrane caps the therapy head** and is pressed **against the patient**.
- A **coupling agent** (gel or oil) is applied between the **rubber membrane** and the **skin to reduce air pockets**.
- **Water inside the therapy head is continuously recirculated and degassed** to prevent bubbles from interfering with **shock wave propagation**.
- **Limitations Compared to Water-Bath Lithotriptors:**
 - **Less effective coupling:**
 - The **rubber membrane** introduces an additional **reflecting interface**.
 - Even with couplants, **tiny air bubbles** between the **skin and membrane** are difficult to eliminate.
 - **In vitro studies show that different couplants significantly affect stone breakage.**
- **Clinical Trade-off:**
 - Despite **reduced performance**, dry therapy heads are **more convenient** for clinical use, making them the **preferred choice in modern lithotriptors**.

The Focal Zone of the Lithotripter

- **Definition and Importance**
 - In lithotripsy, **acoustic energy is focused** to a **small zone** around the **focal point**.
 - The **focal point** is a **geometric location** where the shock waves converge, usually **aligned with the kidney stone for treatment**.
 - In electrohydraulic lithotriptors (EHLs), the **focal point (F2)** is the **second focus of the ellipsoidal reflector**.
 - The **dimensions and pressure characteristics of the focal zone** are key factors **distinguishing different lithotriptors**.
- **IEC Standard Definition of the Focal Zone**
 - Defined as the volume where **peak acoustic pressure is at least half of the maximum peak positive pressure (p^+)**.
 - **Peak positive pressure (p^+)**: The **highest positive pressure** in the waveform (see **Figure 38-1**).
 - **Maximum peak pressure**: The **highest value of p^+ in the lithotripter field**.
 - **Focus location**: Defined as the **point where maximum peak pressure occurs**.
 - The **focal zone shape** is typically an **elongated, elliptical (“cigar-shaped”) volume**.
 - **Important considerations:**
 - The **maximum peak pressure does not always align** with the **stone placement location** recommended by the manufacturer.
 - **Focal zone dimensions change with power settings**.
- **Focal Zone Characteristics in Different Lithotriptors (Figure 38-15)**
 - **Dornier HM3 (Electrohydraulic Lithotripter):**
 - **Maximum peak pressure**: 40 MPa at 20 kV.
 - **Focal zone dimensions**: ~60 mm long, 12 mm diameter.
 - **Storz Modulith (Electromagnetic Lithotripter):**
 - **Maximum peak pressure**: ~100 MPa at energy level 8.
 - **Focal zone dimensions**: ~35 mm long, 4 mm diameter.
 - **Piezoelectric Lithotripter:**
 - **Maximum peak pressure**: ~80 MPa.
 - **Focal zone dimensions**: ~20 mm long, 3 mm diameter.
- **Key Observations**
 - **Significant variation in focal zone dimensions across lithotriptors.**
 - **Lithotriptors with narrower focal zones tend to generate higher peak pressures.**

The Half-Maximum (–6 dB) Focal Zone and Alternative Definitions

- **Standard Definition (IEC Recommendation)**
 - The **half-maximum focal zone** (also called the **−6 dB focal zone**) refers to the **contour where the pressure is 6 decibels lower than the peak pressure**.
 - Example comparisons:
 - **Storz lithotripter** (peak pressure = 110 MPa at energy level 9) → Focal zone defined at 55 MPa.
 - **Dornier HM3** (peak pressure = 40 MPa) → Focal zone defined at 20 MPa.
 - **Problem with this definition:**
 - **Absolute pressure levels differ significantly** between machines.
 - If 55 MPa were used as the threshold for the Dornier HM3, its **focal zone would be nonexistent**.
- **Alternative Focal Zone Definitions**
 - Various alternative definitions have been proposed:
 - **Half the peak negative pressure.**
 - **Half the energy density.**
 - **The surface where peak pressure is 5 MPa.**
 - **Energy passing through a 10 mm diameter volume** (approximate size of a kidney stone).
 - **No universal standard** has been agreed upon due to the **complexity of shock wave fragmentation mechanics**.

Effectiveness of Lithotriptors with Small, Tight Focal Zones

- **Theoretical Advantage**
 - **Electromagnetic (EML) and piezoelectric (PEL) lithotriptors** have **smaller focal spots**, leading to:
 - **More precise targeting** on the stone.
 - **Reduced damage to surrounding tissue**.
- **Experimental and Clinical Findings**
 - **In vitro studies** (with stationary stones) show that:
 - **High-pressure electromagnetic and piezoelectric lithotriptors do not outperform** electrohydraulic lithotriptors (EHLs) in breaking stones.
 - **High peak positive pressure alone does not correlate with better stone fragmentation in clinical settings**.
- **Clinical Issues with Small Focal Zones**
 - **Respiratory-induced stone motion** causes:
 - **Fewer shock waves hitting the stone**.
 - **More energy being deposited into surrounding tissue**.
 - **Lithotriptors with very high peak pressures (>100 MPa)** expose tissue to **excessive acoustic energy**, increasing the risk of **adverse effects** such as:
 - **Subcapsular hematomas**.

Device Equivalency and Equating Lithotripter Performance

- **Lack of Standardized Metrics for Comparison**
 - Currently, **no agreed-upon metrics** exist for comparing the **acoustic output** of different lithotriptors.
 - **No straightforward method** to adjust a lithotripter's settings to **match the performance of another device**.
- **Variability in Lithotripter Characteristics**
 - Although all lithotriptors produce **shock waves with similar waveforms**, key differences include:
 - **Amplitude** (peak positive pressure).
 - **Focal zone size and shape**.
 - **Energy distribution**.
 - **Acoustic field measurements** can yield **very different values** between devices.
- **Example Comparison (Table 38-1)**
 - Comparison of an **electrohydraulic (EHL) vs. electromagnetic (EML) lithotripter**:

- **Energy incident on a 6.5 mm stone:** Similar (~0.484 mJ vs. 0.528 mJ).
 - **Peak positive pressure:** Three times higher in the **electromagnetic lithotripter**.
- **Key takeaway:** Even if one parameter is matched, other properties remain significantly different.
- **Effects of Adjusting Power Settings**
 - Reducing **power** on an **electromagnetic lithotripter** to match the **pressure** of an electrohydraulic lithotripter would cause:
 - **Energy output to drop nearly 10-fold.**
 - **Significant loss in fragmentation efficiency.**
- **Additional Factors Affecting Equivalency**
 - **Number of shock waves:**
 - **Electrohydraulic lithotriptors typically require fewer shock waves** to break stones than electromagnetic lithotriptors.
 - **Shock wave delivery rate:**
 - **The rate of shock wave application** has been reported to **impact fragmentation efficiency.**
- **Conclusion**
 - At present, **no clear way** exists to **adjust power, number of shocks, or delivery rate** to ensure **equivalency** between different lithotriptors.

Waves in Stones

- **Types of Waves in Kidney Stones**
 - Kidney stones are **elastic solids** that support **two types of waves**:
 - Longitudinal (Compression) Waves**
 - Similar to **acoustic waves**.
 - **Molecules oscillate along the direction of propagation** (compression and rarefaction).
 - Travel at **speed c_{Lc_L}** (always faster than transverse waves).
 - Transverse (Shear) Waves**
 - **Oscillations are perpendicular to propagation direction** (similar to a wave on a rope).
 - Unlike longitudinal waves, **do not cause compression and rarefaction**.
 - Travel at **speed c_{Tc_T}** (slower than longitudinal waves).
- **Energy Division at the Stone Boundary**
 - When a **shock wave passes from urine/tissue into a stone**, its energy splits into **longitudinal and transverse waves** inside the stone.
 - **Energy partition depends on:**
 - **Material properties of the stone.**
 - **Angle of incidence** of the wave.
 - **Normal incidence** (perpendicular wave entry):
 - **All energy converts to longitudinal waves.**
 - **No transverse waves generated.**
 - **Oblique incidence** (angled entry):
 - **Less energy enters as longitudinal waves.**
 - **More energy is converted into transverse waves.**
 - **Complex stone shapes** cause a **non-trivial energy distribution** between wave types.
- **Computer Simulation of Shock Wave Interaction (Figure 38-16)**
 - **Simulation models particle motion in an elastic solid.**
 - **Snapshots illustrate stress distribution** inside a cylindrical-shaped stone over time:
 - Initial Frames:**
 - Shock wave enters the stone as **compression waves**.
 - Mid Frames:**

- **Longitudinal waves travel inside the stone.**
- **Acoustic waves outside the stone generate shear waves at the lateral walls.**
- **Shock wave reflects from the rear wall.**

Wave Reflection Effects:

- **Since the fluid impedance is lower than that of the stone, the reflected pressure wave inverts.**
- **The leading compressive wave reflects as a tensile wave (negative pressure reflection coefficient RPR_P due to $Z_2 < Z_1$ $Z_2 < Z_1$).**

Final Frames:

- **Shear and longitudinal waves interact, producing high tensile stresses inside the stone.**
- **These high tensile stresses contribute to stone fragmentation.**

Acoustic Properties of Stones and Shock Wave Fragmentation Mechanisms

Acoustic Properties of Kidney Stones

- The key wave propagation properties of kidney stones include:
 - Density (ρ_0 \rho_0)**
 - Longitudinal sound speed (c_{Lc_L})**
 - Transverse wave velocity (c_{Tc_T})**
- **Variability in Stone Properties (Figure 38-17):**
 - Different stone types (e.g., **uric acid, calcium oxalate monohydrate, cystine stones**) have **widely varying acoustic properties**.
 - Example: **Calcium oxalate monohydrate sound speed** ranges from **3,000 m/s to 4,500 m/s**.
 - Variations may be due to:
 - **Natural differences in stone composition.**
 - **Preparation factors (e.g., hydration levels).**

How Shock Waves Break Stones (Figure 38-18)

- **Several mechanisms have been proposed to explain stone fragmentation by shock waves.**
- Below are some of the **most likely mechanisms**:

Spall Fracture

- **Process:**
 - Shock wave enters the stone and reflects from the rear wall (Figure 38-16).**
 - At the stone/urine interface, the positive pressure pulse inverts.**
 - This results in a large tensile stress near the back wall of the stone.**
 - The incoming negative pressure tail further amplifies the tensile stress.**
 - Since most solids are weaker in tension than compression, this tensile stress causes material failure, leading to fracture.**

Shear Stress Fracture

- **Formation of Shear Stresses:**
 - **Shear waves and compressive waves interact as the shock wave propagates through the stone (Figure 38-19).**
 - Many materials, **especially layered kidney stones**, are **weak in shear**.
 - **Layered structures fracture easily** because:
 - **The bonding strength between layers has low ultimate shear stress.**
 - **The organic binder in kidney stones is softer than the crystalline phase, making it vulnerable to shear forces.**

Shear Stress and Superfocusing in Stone Fragmentation

Shear Stress Fracture

- **Generation of Shear Stresses:**
 - **Shear waves and compressive waves develop as the shock wave propagates through the stone (Figure 38-19).**
 - **Layered kidney stones are particularly weak in shear** because:

- The **bonding strength between layers** has **low ultimate shear stress**.
 - The **organic binder** in kidney stones is **softer than the crystalline phase**.
- As the **shock front moves through the stone**, it induces **high shear stresses at binder/crystal interfaces**, leading to **fracture**.
- **Impact of Shear Waves on Stone Fragmentation:**
 - **Shear waves in the stone can induce tensile stresses** greater than those caused by **spallation**.
 - **Figure 38-16** demonstrates that **shear waves interfere with reflected longitudinal waves**, producing the **largest tensile stresses in cylindrical stones**.

Superfocusing and Stress Amplification

- **What is Superfocusing?**
 - **Superfocusing** occurs when **stress amplification** happens inside the stone due to its **geometry**.
 - The **shock wave reflecting from the distal stone surface** can be focused through:
 - Refraction:**
 - Affected by the **high sound speed and geometry of the stone**.
 - Diffraction:**
 - Occurs at the **stone's corners**, redirecting energy into specific regions.
- **Effects of Superfocusing on Stone Failure:**
 - Reflected waves can **concentrate at caustics** (regions of **high stress** inside the stone).
 - **These high-stress regions lead to material failure**.
 - The **location and intensity of these stresses** depend on the **stone's shape and elastic properties**, such as:
 - **Density** (ρ)
 - **Longitudinal wave speed** (c_{Lc_L})
 - **Shear wave speed** (c_{Tc_T})
- **Conclusion:**
 - **Shear stresses and superfocusing contribute significantly to kidney stone fragmentation**.
 - **Understanding the stone's geometry and wave propagation properties is key to predicting failure locations**.

Squeezing and Cavitation in Stone Fragmentation

Squeezing (Splitting Mechanism)

- **Cause:**
 - Due to the **difference in sound speed** between:
 - **Stone** ($>2,500 > 2,500\text{m/s}$)
 - **Surrounding fluid (urine)** ($\approx 1,500 \approx 1,500\text{m/s}$)
 - This causes the **shock wave inside the stone to travel faster** than the shock wave in the surrounding fluid.
- **Effect on the Stone:**
 - The **fluid shock wave exerts a circumferential force** on the stone (**hoop stress**).
 - This results in **maximum tensile stress** at the **proximal and distal ends** of the stone.
 - Leads to **axial splitting failure** of the stone.
- **Enhancing Squeezing Effect:**
 - The **splitting effect is enhanced** when the **entire stone is within the focal zone**.
 - A **new lithotripter design** has been developed based on this principle and **documented in recent literature**.

Cavitation in Lithotripsy

- **What is Cavitation?**
 - Formation of **small bubbles (cavities)** in the **urine** surrounding the stone due to the **large negative pressure tail of the acoustic pulse**.
- **Microjet Formation and Impact:**
 - When a **cavitation bubble collapses near a solid surface** (e.g., kidney stone):

- A microjet of fluid pierces the bubble and impacts the stone surface at speeds >100 m/s (Figure 38-20).
 - This microjet contributes to stone surface damage.
- **Secondary Shock Wave Emissions:**
 - Bubble collapse also generates secondary shock waves.
 - These waves have amplitudes comparable to the primary shock wave, enhancing fragmentation.
- **Evidence from Experiments:**
 - In vitro studies where cavitation is suppressed show a significant reduction in stone fragmentation.
 - Cavitation is mainly a surface-acting mechanism, primarily affecting:
 - The proximal (shock wave incident) surface of the stone.
- **Cavitation and Spall Fracture:**
 - Cavitation-induced stresses may act via a spallation mechanism.
- **Cluster Cavitation Theory:**
 - Recent research suggests lithotripsy does not produce isolated cavitation bubbles.
 - Instead, it creates bubble clusters (Figure 38-21).
 - The coherent collapse of these clusters enhances cavitation's destructive power.

Fatigue as a Mechanism of Stone Fragmentation

Fatigue Process in Kidney Stones

- **Definition:**
 - Fatigue occurs when progressive crack development leads to stone failure over time.
 - Cracks originate at small imperfections (nucleation sites), which are present in all kidney stones.
- **Formation of Cracks Under Shock Waves:**

Initial imperfections act as stress concentration sites.
Shock waves induce localized stresses far exceeding the average stress.
Repetitive shock waves cause:

 - Imperfections to grow into microcracks.
 - Microcracks to expand into macrocracks.

Tensile or shear stresses eventually enlarge the cracks, leading to failure.
Fatigue is enhanced in regions of high stress coinciding with weaker areas in the stone.
- **Possible Synergistic Effects:**
 - Fatigue may work alongside other mechanisms, such as:
 - Spallation (high tensile stress near the back wall).
 - Shear stress (from longitudinal and transverse wave interactions).
 - Cavitation (localized high-stress regions from bubble collapse).

Evidence Supporting Fatigue as the Main Fragmentation Process

Stone Internal Structure Affects Fragmentation

- Studies show that stone composition and internal microstructure influence how stones break in lithotripsy.

Large Number of Shock Waves Required

- Typically, $>1,000$ shock waves are needed to break stones into small fragments.
- Progressive fracture due to multiple stress cycles is a classic hallmark of fatigue failure.

Uncertainty in the Dominant Fatigue Mechanism

- **Two most commonly cited drivers of fatigue:**

Direct stress mechanisms (spall and shear forces).
Cavitation-induced damage.
Potential combination of both effects.
- **Challenges in Determining the Dominant Mechanism:**
 - Limited material strength data for kidney stones:

- **Compression strength**
- **Fracture toughness**
- **Knoop hardness**
- **Vickers microhardness**
- **Lack of tensile and shear strength data:**
 - **Measuring brittle material properties is technically difficult.**
- **Current mechanical tests are quasi-static** (minutes-long stress applications), whereas **shock waves act in microseconds.**

Conclusion

- **Stone failure in lithotripsy follows a fatigue process, but the dominant mechanism remains unclear.**
- **More data on stone material properties under rapid stress conditions** is needed to fully describe the fracture process.

Mechanisms of Tissue Damage in Shock Wave Lithotripsy (SWL)

Tissue Trauma in Lithotripsy

- **Shock wave lithotripsy (SWL) is known to cause kidney trauma, and in some cases, the injury can be severe.**
- **Clinical implications of SWL-induced damage** are still being investigated (see Chapter 41, “Complications”).
- **Myth: Shock Waves Pass Harmlessly Through the Body**
 - This assumption is **incorrect.**
 - **Most patients receiving at least 2,000 shock waves at midrange power or higher experience some form of tissue trauma.**
 - While **SWL has been highly beneficial**, some patients have suffered **severe or even catastrophic adverse effects.**

Physics of Tissue Damage

- **Acoustic Pulse and Focal Zone Effects**
 - **Lithotriptors generate focused acoustic pulses.**
 - The **acoustic field is broad at the source and narrows at the focus.**
 - The **focal zone has high acoustic pressure and is elongated (~50 mm in length, see Figure 38-15).**
 - Because **kidneys are approximately 50 mm thick**, the **entire kidney is exposed to high-amplitude shock waves.**
 - **Patient motion (respiration, discomfort) can cause the stone to move out of the focal zone, resulting in shock waves impacting tissue instead of the stone.**

Tissue's Resistance to Damage vs. Kidney Stones

- **Advantages of Tissue Properties**
 - **Acoustic impedance of tissue is similar to water**, reducing **shock wave reflection** at tissue-water interfaces.
 - **Tissue does not experience spallation**, as it lacks the **extreme tensile forces** that cause stones to fracture.
 - **Tissue sound speed is nearly constant**, preventing **differential squeezing stresses** that could lead to splitting (unlike stones).

- **Vulnerabilities of Tissue**

Shear Stress

- **Shock wave pressure waves induce shear forces** in tissue.
- Tissue is **not as resistant to shear as compression**, making it susceptible to **mechanical stress-induced injury.**

Cavitation Damage

- The **tensile phase of the shock wave creates cavitation bubbles** in tissue.
- **Bubble collapse generates microjets and secondary shock waves**, which can **cause localized tissue damage.**

Conclusion

- **Tissue is not as vulnerable as kidney stones to certain failure mechanisms (spallation, splitting).**
- **However, it remains susceptible to shear stress and cavitation**, both of which can

contribute to SWL-induced tissue injury.

Mechanical Stress and Cavitation as Mechanisms of Tissue Damage in Lithotripsy

Mechanical Stress in Tissue

- **Compression Due to Positive Pressure**
 - The **positive pressure of a lithotripter pulse** compresses tissue.
 - **Tissue is generally resistant to isotropic compression**, but:
 - The **shock front rise time** is **~70 ns** in tissue.
 - This corresponds to a **spatial scale of ~100 μm** .
 - Tissue structures between **10 μm to 1 mm** experience **large stress variations** as the shock wave passes.
- **Shear Forces from Shock Waves**
 - The **short rise time of the shock wave** results in **non-uniform tissue straining**, leading to **shear forces**.
 - **Tissue structures are highly sensitive to shear stress**.
 - **Shock wave-induced tissue distortion** could produce **enough shear to cause mechanical damage**.

Shear Stress Due to Tissue Inhomogeneities

- **Tissue is an inhomogeneous medium** with varying sound speeds at different length scales.
- **Sound speed variations (~mm scale) affect wave focusing**:
 - **Parts of the shock wave traveling through higher sound-speed regions** arrive **sooner**.
 - **Parts traveling through lower sound-speed regions** lag behind.
 - This **distortion of the wavefront** leads to **localized shear stresses in tissue**.
 - **Shear stress could be strong enough to cause mechanical damage**.

Cavitation as a Dominant Cause of Tissue Damage

- **Cavitation in Tissue During Lithotripsy**
 - **Cavitation bubbles form in tissue** due to the **tensile phase of the shock wave**.
 - **Passive cavitation detection (PCD) studies in humans and pigs** have identified **cavitation signatures in**:
 - **Perirenal fat**.
 - **The collecting system**.
 - **Renal parenchyma**.
 - **Subcapsular hematomas**.
- **Cavitation's Biological Effects**
 - Cavitation has been well-documented to cause **biological damage in vitro** (see **Carstensen EL, Miller DL, Dalecki D, Delius M, Williams JC**).
 - **Experimental evidence**:
 - **Reducing cavitation significantly reduces tissue damage** in both **in vitro** cells and **in vivo** tissue.
- **Cavitation in Blood Vessels vs. Tissue**
 - **Cavitation is more damaging in blood vessels than in surrounding tissue**:
 - **Tissue-constrained bubbles** cannot expand or collapse violently.
 - **Bubbles in blood vessels** experience **unrestricted growth and collapse**, making them **more destructive**.
 - This suggests that **vascular injury is a key concern in lithotripsy-induced tissue damage**.

Cavitation-Induced Tissue Damage Mechanisms

1. Collapsing Bubbles and Microjets

- **Mechanism**:
 - **Asymmetrical bubble collapse** generates **high-velocity microjets** of fluid.
 - These microjets are **strong enough to pit foils or etch metal**.
 - **Blood vessels, particularly capillaries**, could be **punctured** by these jets.

- **Weakness of the Theory:**
 - **Most injured blood vessels during lithotripsy are too small** to allow for **full bubble expansion and collapse**.
 - This suggests that **collapse-induced damage may not fully explain the vascular injuries seen in SWL**.

2. Bubble Expansion and Vessel Rupture

- **Mechanism:**
 - As the **shock wave's negative pressure phase passes through a blood vessel**, bubbles experience **rapid expansion** (see **Figure 38-10**).
 - **Explosive bubble growth exerts outward pressure**, potentially **rupturing vessel walls**.
- **Supporting Evidence:**
 - **Capillary rupture is often observed first**, indicating **small vessels experience the greatest stress**.
 - **In vitro experiments using capillary phantoms** support the **explosive bubble hypothesis**.
- **Potential for Additional Tissue Damage:**
 - If **bubble expansion can rupture blood vessels**, it may also **tear apart other soft tissue structures nearby**.

3. Hematoma Formation and Persistent Cavitation

- **Once blood vessels rupture, blood pools** in areas like hematomas, creating:
 - **Larger fluid-filled spaces** for cavitation bubbles to form.
 - **Pre-existing bubbles** in these pools that **act as nuclei for further cavitation events**.
 - **Reduced blood flow**, preventing cavitation bubbles from being **swept away**, prolonging their effects.
- **Observed Effects in SWL:**
 - **Intense passive cavitation detection (PCD) signals** recorded from hematomas.
 - **B-scan ultrasound shows increased echogenicity**, indicating **cavitation activity within hematomas**.
 - **Violent cavitation in pooled blood may lead to additional cell disruption** in the affected region.

Current Research and Unresolved Questions

- **Understanding tissue response to shock waves is still incomplete:**
 - **Tissue mechanical properties at high strain rates** (relevant to SWL) are **poorly understood**.
 - **Damage criteria are not well-defined**, making it **difficult to quantify tissue injury risk**.
- **Few experimental models exist to validate different hypotheses** about cavitation-induced damage.
- **General consensus:**
 - **Cavitation is likely the primary mechanism of tissue injury in SWL**.
 - **Further research is needed** to confirm the dominant physical processes involved.

The Evolution of the Lithotripter

How Lithotriptors Have Changed Over the Years

- **Lithotripsy has been in clinical practice for over 25 years**, with several design improvements, but **no fundamental changes in acoustics**.
- **Key advancements in modern lithotriptors:**
 - **Compact and modular designs**.
 - **Dry shock heads** instead of water-bath immersion.
 - **Improved imaging for stone targeting**.
- **Despite these changes, the acoustic signature of lithotriptors remains largely the same.**
 - The **focal waveform of the Dornier HM3 is nearly identical** to that of modern lithotriptors.

Trade-Offs in Modern Lithotripsy

- **Success rates in lithotripsy have declined** with the introduction of newer devices.
 - **Machines with tight focal zones and extreme peak positive pressures** have:
 - **Reduced stone fragmentation efficiency.**
 - **Increased collateral tissue damage.**
 - **Progress has been made in convenience, but at a cost to treatment effectiveness.**
- Early Lithotriptors: Electrohydraulic Devices**
- The **first-generation lithotriptors** used **electrohydraulic shock wave generation**:
 - **Shock waves were produced by an underwater spark discharge.**
 - **Patients were fully immersed in a water bath** for optimal shock wave coupling.

The Dornier HM3: A Landmark Lithotripter

- The **Dornier HM3** was one of the most successful **first-generation** lithotriptors.
- **Key characteristics of the HM3**:
 - **Moderate peak positive pressures (~40 MPa).**
 - **Generous focal zone (~15 mm by 60 mm).**
- **Impact of the Dornier HM3**:
 - **Widely adopted and highly successful.**
 - **Played a crucial role in the early success and widespread acceptance of lithotripsy.**
 - **Some centers still use the HM3 today** due to its **proven effectiveness.**

Challenges and Advancements in Lithotripter Design

Perceived Drawbacks of the Dornier HM3

- Despite its effectiveness, the **Dornier HM3** had several limitations:
 - Open water bath** was required for shock wave coupling.
 - Treatment was painful**, requiring **sedation or anesthesia.**
 - Large, stationary equipment** required a **dedicated water treatment facility.**
- **Physicians and Patients Desired Improvements**:
 - **A painless, convenient, walk-in/walk-out therapy** with **minimal or no anesthesia.**
 - Lithotripter manufacturers **modified designs to address these concerns.**

Modifications and Improvements in Modern Lithotriptors

Elimination of the Open Water Bath

- **Shock head enclosed with a rubber membrane** for body coupling.
- **Advantages**:
 - **Easier access to patients** for medical staff.
 - **Lithotriptors no longer required a dedicated facility.**
 - **Enabled modular, portable lithotripter designs** (used in mobile lithotripsy units).
- **Disadvantages**:
 - **Acoustic coupling with a water-tissue interface is superior.**
 - **Membrane coupling is less effective**, leading to **energy loss.**

Replacement of Electrohydraulic Electrodes

- **Electrohydraulic lithotriptors** required **frequent electrode replacements**, increasing cost and treatment time.
- **Electromagnetic and piezoelectric lithotriptors** **do not use electrodes**, reducing maintenance.
- **Issues with electrohydraulic electrode wear**:
 - **Spark gap widens with use**, leading to **greater variability in arc discharge paths.**
 - **Manufacturers improved electrode lifespan** using **electrolyte-filled encapsulated electrodes.**

The Challenge of Anesthesia-Free Lithotripsy

Discomfort from Shock Wave Treatment

- **Main cause of pain**: **Cutaneous pain at the shock wave entry point.**
- **Attempts to reduce pain** included **widening the shock source aperture** to spread energy over a broader area.

Trade-Offs of a Wider Aperture

- Broader acoustic field along the shock wave axis.
- Narrower focal zone (~5 mm or less).
- Some modern lithotriptors generate extreme peak positive pressures (>100 MPa).
- Problems with a tight focal zone:
 - Difficult to keep perfectly aligned with the stone due to respiratory motion.
 - More challenging to hit the stone compared to a broad focal zone.

Current Reality of Lithotripsy Comfort

- Regardless of design improvements, lithotripsy remains uncomfortable.
- If the patient is **not sedated**, they will **move during treatment**, reducing effectiveness.
- Attempts to create a completely anesthesia-free lithotripter have not yet been successful.

Future Directions in Lithotripter Design

1. Wide-Focus, Low-Pressure Lithotriptors

- New approach in response to high-pressure, tight-focal-zone machines.
- Xi Xin-Eisenmenger Lithotripter:
 - Largest focal zone in clinical use (18 × 180 mm).
 - Low acoustic pressure (10–25 MPa).
 - Designed to test the **circumferential squeezing hypothesis** for enhanced stone breakage.
 - Early trial results:
 - High stone-free rate (86%).
 - Can be used without anesthesia.

2. Cavitation Control Strategies

- Manipulating cavitation bubbles with a second pulse:
 - Cavitation bubble cycle lasts ~300 μs in free field, ~600 μs at a stone surface.
 - Effect of timing on bubble collapse:
 - Early second pulse: Stops further expansion → **minimal damage**.
 - Late second pulse: Accelerates collapse → **enhanced damage**.
- Dual Pulse Lithotripsy:
 - Developed by Bailey using **twin shock sources facing each other**.
 - Other researchers have refined this by:
 - Firing multiple pulses along the same axis.
 - Using dual treatment heads at an angle to accommodate patient anatomy.
 - Potential benefits:
 - Tailored acoustic forces for improved stone breakage.
 - Reduced collateral tissue damage.
 - Still under development and testing—no definitive efficacy and safety data yet.

3. Real-Time Stone Tracking for Targeted Shock Wave Delivery

- Goal: Reduce unnecessary shock wave impact on tissue.
- Proposed tracking methods:
 - Monitor stone location and fire only when the stone is at the focus.
 - Acoustic time-reversal technology:
 - Dynamically adjusts lithotripter focus to track the stone as it moves.
- Potential impact:
 - Fewer shock waves needed to break stones.
 - Lower risk of tissue injury.
- Current limitation:
 - No clinical devices yet employ real-time tracking.

Final Thoughts

- Lithotripsy technology is evolving toward improving efficiency and reducing tissue damage.

- **Wide-focus, low-pressure lithotriptors, dual-pulse strategies, and real-time tracking** hold promise.
- **Further research and clinical validation are needed** before these innovations become standard practice.